

16S Pipeline: Software Requirements Specification

Author: Sarah Burns

Date: 14/10/2019

1. Sign offs

I agree that this document represents our best understanding of the requirements for this project today and the system described will satisfy our needs.

Meera Chand

Name

Signature

Date

Penny Cliff

Name

Signature

Date

Rossella Baldan

Name

Signature

Date

2. Introduction

2.1. Purpose

To develop an analytical pipeline to detect bacterial species from 16S nanopore sequencing data.

2.2. Scope

The pipeline will be developed within the Centre for Clinical Infection and Diagnostics Research (CIDR) with the eventual aim of inclusion within a UKAS-accredited clinical service. The analysis will be automated and performed in real-time, and be managed by a bioinformatician in addition to selected personnel.

2.3. Definitions, acronyms and abbreviations

ATCC	American Type Culture Collection (reference standards)
DNAexus	Cloud-based data analysis and management platform
WDL	Workflow development language; a framework for writing pipelines.
Nextflow	A framework for writing pipelines.

3. Overall description

3.1. Product perspective

This will form part of a new biological diagnostic assay with improved speed and accuracy compared to current gold-standard methods. The pipeline will evolve over the course of a pilot study prior to implementation as a routine clinical service.

3.2. Product features

The raw output of the nanopore sequencing is FAST5 files which are produced and base-called in real-time, producing FASTQ files. The pipeline will monitor the generation of the FASTQ files and, at select time intervals, merge these into a single FASTQ file, classify the reads against a reference database and produce a clinically useful report alongside quality control information (e.g. read length, yield, base-call quality). See Figure 1 in Supporting information.

3.3. Operating environment

- Pipeline framework: DNAnexus-compatible (e.g. WDL or Nextflow)
- Operating system: Linux

3.4. Assumptions and dependencies

- A maximum of five samples will be sequenced in parallel.
- All latest versions of the necessary software will be installed in the run environment.
- Input data will be in the correct format.
- There will be a continuous internet connection.

4. System features

4.1. Functional requirements

<u>Requirement</u>	<u>Priority</u>
Database to include sequences for all relevant clinical species (see Table 1 in Supporting information).	<i>Essential</i>
Classifier can detect to species level.	<i>Essential</i>
Pipeline output to report hits in >1% classified reads (subject to testing).	<i>Essential</i>
Quality report to be produced for every sample.	<i>Essential</i>
Pipeline executable both within DNAnexus and on local workstation.	<i>Desirable</i>
Data to be stored within DNAnexus.	<i>Desirable</i>

4.2. Non-functional requirements

<u>Requirement</u>	<u>Priority</u>
Capacity Able to process five samples in parallel.	<i>Essential</i>
Performance Complete within 30 minutes (larger files >1GB permitted to exceed this providing the duration remains within an acceptable timeframe for rapid clinical service).	<i>Desirable</i>
Accessibility Continuous service.	<i>Desirable</i>

Maintainability and extensibility	Code will be modular, easy to understand and written in line with conventions, with detailed commenting and effective error handling.	<i>Desirable</i>
-----------------------------------	---------------------------------------------------------------------------------------------------------------------------------------	------------------

4.3. Software releases

Initial release (v1.0.0) containing all essential requirements as a minimum will be available prior to the start of the pilot study.

4.4. Validation

- 167 ATCC isolates to be used as “truth” data to validate classifier.
- Clinical samples containing known species to be used to determine sensitivity and specificity of test.

4.5. Documentation

- Standard operating procedure
- Testing specification
- Validation specification

5. Supporting information

Figure 1 – General pipeline workflow




Table 1 – Clinical species desirable for pipeline to detect

<i>Abiotrophia defectiva</i>	<i>Coxiella burnetii</i>	<i>Peptoniphilus harei</i>
<i>Achromobacter xylosoxidans</i>	<i>Cutibacterium acnes</i>	<i>Peptostreptococcus anaerobius</i>
<i>Acinetobacter baumannii</i>	<i>Cutibacterium avidum</i>	<i>Plesiomonas shigelloides</i>
<i>Acinetobacter baylyi</i>	<i>Delftia acidovorans</i>	<i>Pluralibacter gergoviae</i>
<i>Acinetobacter haemolyticus</i>	<i>Dermabacter hominis</i>	<i>Prevotella denticola</i>
<i>Acinetobacter johnsonii</i>	<i>Eikenella corrodens</i>	<i>Proteus mirabilis</i>
<i>Acinetobacter junii</i>	<i>Elizabethkingia miricola</i>	<i>Proteus vulgaris</i>
<i>Acinetobacter lwoffii</i>	<i>Enterobacter asburiae</i>	<i>Providencia rettgeri</i>
<i>Acinetobacter nosocomialis</i>	<i>Enterobacter cloacae</i>	<i>Providencia stuartii</i>
<i>Acinetobacter pittii</i>	<i>Enterobacter kobei</i>	<i>Pseudomonas aeruginosa</i>
<i>Acinetobacter radioresistens</i>	<i>Enterobacter ludwigii</i>	<i>Pseudomonas koreensis</i>
<i>Acinetobacter tjernbergiae</i>	<i>Enterococcus avium</i>	<i>Pseudomonas luteola</i>
<i>Acinetobacter ursingii</i>	<i>Enterococcus casseliflavus</i>	<i>Pseudomonas mosseli</i>
<i>Actinomyces naeslundii</i>	<i>Enterococcus durans</i>	<i>Pseudomonas oryzihabitans</i>
<i>Actinomyces neuui</i>	<i>Enterococcus faecalis</i>	<i>Pseudomonas putida</i>
<i>Actinomyces odontolyticus</i>	<i>Enterococcus faecium</i>	<i>Pseudomonas stutzeri</i>
<i>Actinomyces radingae</i>	<i>Enterococcus raffinosus</i>	<i>Pseudomonas taiwanensis</i>
<i>Actinomyces turicensis</i>	<i>Erysipelotrix rhusiopathiae</i>	<i>Raoultella ornithinolytica</i>
<i>Aerococcus urinae</i>	<i>Escherichia coli</i>	<i>Raoultella planticola</i>
<i>Aeromonas caviae</i>	<i>Finnegoldia magna</i>	<i>Rhodococcus hoagii</i>
<i>Aeromonas media</i>	<i>Fusobacterium necrophorum</i>	<i>Salmonella enterica</i>
<i>Aeromonas veronii</i>	<i>Fusobacterium nucleatum</i>	<i>Serratia liquefaciens</i>
<i>Aggregatibacter actinomycetemcomitans</i>	<i>Gemella haemolysans</i>	<i>Serratia marcescens</i>
<i>Aggregatibacter aphrophilus</i>	<i>Gemella morbillorum</i>	<i>Serratia ureilytica</i>
<i>Aggregatibacter segnis</i>	<i>Gemella sanguinis</i>	<i>Staphylococcus aureus</i>
<i>Agrobacterium tumefaciens</i>	<i>Granulicatella adiacens</i>	<i>Staphylococcus capitis</i>
<i>Alcaligenes faecalis</i>	<i>Haemophilus haemolyticus</i>	<i>Staphylococcus caprae</i>
<i>Anaerococcus lactolyticus</i>	<i>Haemophilus influenzae</i>	<i>Staphylococcus epidermidis</i>
<i>Anaerococcus octavius</i>	<i>Haemophilus parainfluenzae</i>	<i>Staphylococcus haemolyticus</i>
<i>Arcanobacterium haemolyticum</i>	<i>Hafnia alvei</i>	<i>Staphylococcus hominis</i>
<i>Bacillus cereus</i>	<i>Helcococcus kunzii</i>	<i>Staphylococcus lugdunensis</i>
<i>Bacillus licheniformis</i>	<i>Kingella kingae</i>	<i>Staphylococcus pasteurii</i>
<i>Bacillus subtilis</i>	<i>Klebsiella aerogenes</i>	<i>Staphylococcus pettenkoferi</i>
<i>Bacteroides caccae</i>	<i>Klebsiella oxytoca</i>	<i>Staphylococcus pseudointermedius</i>
<i>Bacteroides fragilis</i>	<i>Klebsiella pneumoniae</i>	<i>Staphylococcus saprophyticus</i>
<i>Bacteroides thetaotaomicron</i>	<i>Klebsiella variicola</i>	<i>Staphylococcus schleiferi</i>
<i>Bartonella henselae</i>	<i>Lactobacillus gasseri</i>	<i>Staphylococcus simulans</i>
<i>Bartonella quintana</i>	<i>Lactobacillus rhamnosus</i>	<i>Staphylococcus warneri</i>
<i>Bordetella pertussis</i>	<i>Lautropia mirabilis</i>	<i>Stenotrophomonas maltophilia</i>
<i>Brevibacterium casei</i>	<i>Legionella longbeachae</i>	<i>Streptococcus agalactiae</i>
<i>Brucella abortus</i>	<i>Legionella pneumophila</i>	<i>Streptococcus anginosus</i>
<i>Burkholderia cenocepacia</i>	<i>Listeria monocytogenes</i>	<i>Streptococcus constellatus</i>
<i>Burkholderia cepacia</i>	<i>Massilia timonae</i>	<i>Streptococcus cristatus</i>
<i>Burkholderia gladioli</i>	<i>Micrococcus luteus</i>	<i>Streptococcus dysgalactiae</i>
<i>Cardiobacterium hominis</i>	<i>Moraxella catarrhalis</i>	<i>Streptococcus gallolyticus</i>
<i>Chlamydomphila pneumoniae</i>	<i>Moraxella nonliquefaciens</i>	<i>Streptococcus gordonii</i>
<i>Citrobacter amalonaticus</i>	<i>Moraxella osloensis</i>	<i>Streptococcus intermedius</i>
<i>Citrobacter braakii</i>	<i>Morganella morganii</i>	<i>Streptococcus lutetiensis</i>
<i>Citrobacter freundii</i>	<i>Mycobacterium abscessus</i>	<i>Streptococcus mitis</i>
<i>Citrobacter koseri</i>	<i>Mycobacterium avium</i>	<i>Streptococcus mutans</i>
<i>Clostridium paraputrificum</i>	<i>Mycobacterium chimaerae</i>	<i>Streptococcus oralis</i>
<i>Clostridium perfringens</i>	<i>Mycobacterium fortuitum</i>	<i>Streptococcus parasanguinis</i>
<i>Comamonas testosteroni</i>	<i>Mycobacterium intracellulare</i>	<i>Streptococcus pastorianus</i>
<i>Comamonas kerstersii</i>	<i>Mycoplasma lipophilum</i>	<i>Streptococcus pneumoniae</i>
<i>Corynebacterium amycolatum</i>	<i>Mycoplasma pneumoniae</i>	<i>Streptococcus pyogenes</i>
<i>Corynebacterium aurimucosum</i>	<i>Neisseria gonorrhoea</i>	<i>Streptococcus salivarius</i>
<i>Corynebacterium confusum</i>	<i>Neisseria meningitidis</i>	<i>Streptococcus sanguinis</i>
<i>Corynebacterium diphtheriae</i>	<i>Neisseria weaveri</i>	<i>Sutterella wadsworthensis</i>
<i>Corynebacterium haemolyticum</i>	<i>Ochrobactrum anthropi</i>	<i>Tropheryma whipplei</i>
<i>Corynebacterium jeikeium</i>	<i>Oligella urethralis</i>	<i>Trueperella bernardiae</i>
<i>Corynebacterium minutissimum</i>	<i>Paenibacillus illinoisensis</i>	<i>Veillonella atypica</i>
<i>Corynebacterium propinquum</i>	<i>Paenibacillus pabuli</i>	<i>Veillonella parvula</i>
<i>Corynebacterium striatum</i>	<i>Pantoea agglomerans</i>	<i>Vibrio fluvialis</i>
<i>Corynebacterium tuberculoostearicum</i>	<i>Parvimonas micra</i>	<i>Vibrio parahaemolyticus</i>
<i>Corynebacterium ulcerans</i>	<i>Pasteurella multocida</i>	